


Objective: Determine the concentration of Nickel and Cobalt in an unknown solution using ion-exchange separation: Nickel will be complexed to form a red-brown complex, and cobalt will be converted to blue $Co(CNS)_4^{2-}$; both will be measured spectrophotometrically


theory

Ion exchange equilibria: treat with law of mass exchange.



for which
$$K = \frac{[(res^+) x A^{x-}] [Cl^- (aq)]^x}{[(res^+) Cl^-]^x [A^{x-} (aq)]}$$

For $CoCl_4^{2-}$ ions and basic ion exchange resin $\phi - CH_2 N^+(CH_3)_3 Cl^-$



with large K such that $CoCl_4^{2-}$ is retained on the resin @ low pH.

ions like Ni^{2+} that do not complex with Cl^- , K is small and Ni^{2+} is not retained (thus separating the Ni from the Co)

* Ni red complex: oxidize Ni w/ bromine in ammoniacal solution + treat with dimethylglyoxime (measure @ 450 nm)


* Co blue complex: convert metal to $Co(CNS)_4^{2-}$ in acetone/ethanol solution to prevent complex dissociation in H_2O (measure @ 625 nm)

Beer - Lambert law

$$A = \epsilon c l$$

↑ ↑
y x

$$m = \epsilon l$$





The introduction/goals/theory section should be ~1 page for an in-person experiment and ~3 pages for an online lab

Exp. No. 12	Experiment/Subject Ion-exchange separation	Date 3/19/20
Name	Lab Partner	Locker/Desk No. Course & Section No. 118

Your flow chart or abbreviated steps only need to be as detailed as you need to perform the experiment by using them without having to consult the experimental writeup.

Procedure



Signature sign	Date date	Witness/TA	Date
-------------------	--------------	------------	------

Exp. No. 12	Experiment/Subject Ion-Exchange Separation	Date 3/19/20
Name	Lab Partner	Locker/ Desk No.
		Course & Section No. 118

c)

Unknown: Ni + Co in 9M HCl

add 900 0.2 mL with 250 μ l pipettor

- but tip as close to resin bed as possible
- add solution as slowly as possible

Clean 250 mL volumetric flask

solution seeps into resin
blue (cobalt) band

- will spread during washing but will not come off column

4 1-mL portions 9M HCl (allow each to settle before adding next portion)

end of last wash

do not scratch out

catch ~~some~~ drop in beaker

add 1 drop concentrated NH_3
+ 1 drop dimethyl glyoxime
- no red color should form

4-5 1 mL portions 1M HCl


band should move down + into flask
(color blue first, then pink)

Ni flask (set aside)


50 mL vol. flask

Signature sign	Date date	Witness/TA	Date
-------------------	--------------	------------	------


Exp. No. 12	Experiment/Subject 16h - Exchange separation	Date 3/19/20	
Name	Lab Partner	Locker/ Desk No.	Course & Section No.



* to each standard flask




- measure absorbance at each solution, using blank as reference (do not wait more than 30 min)
- measure @ 450 nm



- 88 mL 95% ethanol
- 50 mL dimethylglyoxime
- dilute to mark w/ H_2O + measure absorbance in cuvette @ 450 nm


Signature sign	Date date	Witness/TA	Date
-------------------	--------------	------------	------

Exp. No. 12	Experiment/Subject ion-exchange separation	Date 3/19/2020
Nam	Lab Partner	Locker/ Desk No.
		Course & Section No. 118



F) Calculations

plot graph of A vs concentration for Ni , Co standard solutions



using least squares method

$$A = \frac{Ebc}{m \times x}$$

→ use to calculate $[\text{Ni}]$ and $[\text{Co}]$ from unknown

calculate uncertainty:

$$u_x = \frac{S_y}{|m|} \sqrt{\frac{1}{k} + \frac{1}{n} + \frac{(y - \bar{y})^2}{m^2 \sum (x_i - \bar{x})^2}}$$

where $n, S_y, m, x_i, \bar{x}, \bar{y}$ refer to standard
 k, y refer to unknowns

calculate concentrations in mg/ml

Signature sign	Date date	Witness/TA	Date
-------------------	--------------	------------	------