

PEXMUX

- PEXMUX 8/16 Control Unit
User Manual
Version 002



Copyright © by Bruker Corporation

All rights reserved. No part of this publication may be reproduced, stored in a retrieval system, or transmitted, in any form, or by any means without the prior consent of the publisher. Product names used are trademarks or registered trademarks of their respective holders.

This manual was written by

Willy Uhrig and Thomas Eckert

© August 16, 2013: Bruker Corporation

Rheinstetten, Germany

P/N: Z31929

DWG-Nr.: Z4D11282A

For further technical assistance on the PEXMUX unit, please do not hesitate to contact your nearest BRUKER dealer or contact us directly at:

BRUKER Corporation
am Silberstreifen
D-76287 Rheinstetten
Germany

Phone: + 49 721 5161 0

FAX: + 49 721 5171 01

E-mail: nmr-support@bruker.de

Internet: www.bruker.com

1	About	5
1.1	This Manual	5
1.2	Policy Statement	5
1.3	Symbols and Conventions.....	5
2	Introduction.....	7
2.1	Concept	7
2.2	Limitation of Liability	7
2.3	Before You Begin	7
2.4	Spare Parts	8
2.5	Warranty Terms.....	8
2.6	The Bruker Service.....	8
3	Safety.....	9
3.1	Intended Use	9
3.2	General Safety Instructions.....	9
4	Technical Data	11
4.1	General Information.....	11
4.2	Utility Requirements	11
4.3	Operating Conditions.....	12
4.3.1	Operating Environment	12
4.4	Rating Plates	13
5	Design and Function	15
5.1	Product and Manufacturer.....	15
5.2	General Description.....	15
5.3	Features	16
5.4	Functional Description.....	17
5.5	The PEXMUX Board H12583.....	20
5.6	JTAG Structure.....	20
5.7	Power Requirements.....	22
5.8	Pin Allocation for the LVDS Connectors.....	23
5.8.1	CHASSIS.....	23
6	Transport, Packaging and Storage	27
6.1	Inspection at Delivery	27
6.2	Packaging.....	27
6.2.1	Handling Packaging Materials.....	27
6.3	Storage.....	28
6.4	Transport to Manufacturer.....	28

7	Installation and Initial Commissioning	29
7.1	Safety	29
7.2	Installation	30
7.2.1	External Connections	31
8	Operating Instructions	33
8.1	Operator Protection	33
8.2	Unit Protection	34
9	Maintenance and Cleaning	35
9.1	Safety	35
9.1.1	Environmental Protection	36
9.2	Maintenance	36
9.3	Cleaning	36
9.4	Service Requests	36
10	Dismantling and Disposal	37
11	Contact	39
Appendix A		41
A.1	Warning Signs	41
A.2	Figures.....	43
A.3	Tables.....	45

1 About

1.1 This Manual

This manual is intended to be a reference guide for operators and service technicians. It provides detailed information about the user level maintenance and service and overall use of the Bruker device.

This manual is an integral part of the device, and must be kept in close proximity to the device where it is permanently accessible to personnel.

Before starting any work, the personnel must have read the manual thoroughly and understood its contents. Compliance with all specified safety instructions and operating instructions is vital to ensure safe operation.

In addition, local accident prevention regulations and general safety instructions must be observed for the operational area of the device.

Illustrations in this manual are intended to facilitate basic understanding, and may differ from the actual design.

In addition to the user manual, instructions concerning local labor protection laws, operator regulations, and tools and supplies must be available and adhered to.

1.2 Policy Statement

It is the policy of Bruker to improve products as new techniques and components become available. Bruker reserves the right to change specifications at any time.

Every effort has been made to avoid errors in text and figure presentation in this publication. In order to produce useful and appropriate documentation, we welcome your comments on this publication. Support engineers are advised to regularly check with Bruker for updated information.

Bruker is committed to providing customers with inventive, high quality products and services that are environmentally sound.

1.3 Symbols and Conventions

Safety instructions in this manual are marked with symbols. The safety instructions are introduced using indicative words which express the extent of the hazard.

In order to avoid accidents, personal injury or damage to property, always observe safety instructions and proceed with care.



DANGER

This combination of symbol and signal word indicates an immediately hazardous situation which could result in death or serious injury unless avoided.



WARNING

This combination of symbol and signal word indicates a potentially hazardous situation which could result in death or serious injury unless avoided.



CAUTION

This combination of symbol and signal word indicates a possibly hazardous situation which could result in minor or slight injury unless avoided.

SAFETY INSTRUCTIONS

This combination of color and signal words are used for control flow and shutdowns in the event of an error or emergency.

NOTICE

This combination of color and signal word indicates a possibly hazardous situation which could result in damage to property or the environment unless avoided.



This symbol highlights useful tips and recommendations as well as information designed to ensure efficient and smooth operation.

2 Introduction

2.1 Concept

The Parallel Excitation Multiplexing (PEXMUX) Unit is used to split the available IPSO TX-Controller LVDS channels in **AVANCE III** spectrometers.

This manual is intended to be used by trained device users. It contains information about the device including operation, safety, and maintenance.

2.2 Limitation of Liability

All specifications and instructions in this manual have been compiled taking account of applicable standards and regulations, the current state of technology and the experience and insights we have gained over the years.

The manufacturer accepts no liability for damage due to:

- Failure to observe this manual.
- Improper use.
- Deployment of untrained personnel.
- Unauthorized modifications.
- Technical modifications.
- Use of unauthorized spare parts.

The actual scope of supply may differ from the explanations and depictions in this manual in the case of special designs, take-up of additional ordering options, or as a result of the latest technical modifications.

The undertakings agreed in the supply contract as well as the manufacturer's Terms and Conditions and Terms of Delivery and the legal regulations applicable at the time of conclusion of the contract shall apply.

2.3 Before You Begin

This user manual contains information and safety information that are necessary for the safe operation of the device.

Any user maintenance and repairs are to be accomplished using the information in this manual.

Consider all safety references!

Information for ordering spare parts is available in the spare parts section for from the Bruker Service Center (see contacts).

2.4 Spare Parts

Spare parts may be obtained from authorized dealers or directly from the manufacturer. See "[Contact](#)" on page 39 for the address.

2.5 Warranty Terms

The warranty terms are included in the manufacturer's Terms and Conditions.

2.6 The Bruker Service

Our customer service division is available to provide technical information. See "[Contact](#)" on page 39 for contact details.

In addition, our employees are always interested in acquiring new information and experience gained from practical application; such information and experience may help improve our products.

3 Safety

This chapter provides an overview of all the main safety aspects involved in ensuring optimal personnel protection and safe and smooth operation.

Non-compliance with the action guidelines and safety instructions contained in this manual may result in serious hazards. Before starting any work, personnel must read this chapter thoroughly and understand its contents.

3.1 Intended Use

The unit has been designed and constructed solely for the intended use described here.

The PEXMUX units must only be used for the limited purpose of splitting the available IPSO TX-controller LVDS channels in Bruker AVANCE III spectrometers.

Intended use also includes compliance with all specifications in this manual.

Any use which exceeds or differs from the intended use shall be considered improper use.

No claims of any kind for damage will be entertained if such claims result from improper use.

3.2 General Safety Instructions

The PEXMUX can be damaged by inappropriate usage. In this case the equipment must not be used until it has been checked service personnel.

The user should re-check the equipment at regular intervals for any damage or wear and is expected to inform the service immediately of any abnormality.

Do not operate the equipment in the presence of flammable gases or fumes.

WARNING

Material or personnel damage from worn or damaged parts.

Damaged or worn parts may result if one or the following circumstances occurs:

- ▶ The power cord, power plug or power supply are cracked, brittle or damaged.
- ▶ Signs of excessive heat appear.
- ▶ There is evidence or suspicion that a liquid has intruded into any enclosure.
- ▶ The power cord or the power supply have been in contact with any liquid.
- ▶ The PEXMUX has been dropped or damaged in any way.
- ▶ The equipment does not work correctly

In the unlikely event one of the above occurs:

- ▶ Stop using the equipment.
- ▶ Disconnect the power supply.
- ▶ Inform service and ask for instructions.



i Only trained service personnel are allowed to mount, retrofit, repair, adjust and dismantle the PEXMUX. Do not try to service the equipment unless you are specifically asked to do so and are given instructions by the service staff. In case of questions or problems, please contact your nearest Bruker office or representative.

i Before maintenance, repair or shipment, the PEXMUX Unit must be completely disconnected from the power supply and dismantled from its rack.

4 Technical Data

4.1 General Information

Specification	Value	Unit
Weight	2.5	kg
Length	46.2	cm
Width	29.0	cm
Height	4.3	cm

4.2 Utility Requirements

Ingress Protection Class:	IP 20, secured against touching dangerous points inside, not secured against ingress of water.
Over-voltage Protection Category:	CAT II, according to EN61010–1 or IEC60664–1 safe against over-voltage by switching, not safe against lightning. Surge immunity of PSU is level 3, according to IEC61000–4–5.
IEC Protection Class:	Class I, according the IEC 61140, connected to mains by phase, neutral and protective earth.
Input Voltage Range:	208 V to 230 V.
Frequency:	50/60 Hz.
Input Power:	30 VA.
Fuse On Mains Input:	250 VAC 2 A.
Connection:	Socket outlet with phase, neutral and protective-earth according to VDE 0620–1.

4.3 Operating Conditions

4.3.1 Operating Environment

Permissible ambient temperature:	5°C to 40°C.
Permissible altitude:	Up to 2000 meters above sea level.
Relative humidity:	Maximum of 80% for temperatures up to 31°C, linearly decreasing to 50% at a temperature of 40°C.
Permissible storage temperature:	5°C to 40°C.
Pollution degree:	The acceptable pollution degree is 2. According to EN 61010-1 or IEC 60664-1. Any pollution is non-conductive, except condensed moisture.

4.4 Rating Plates



The rating plate is located at the power input and includes the following information:

- Manufacturer
- Type
- Voltage
- Frequency
- Apparent power consumption, maximum
- PN: Part Number
- SN: Serial Number
- Va: Variant
- ECL: Engineering Change Level



5 Design and Function

5.1 Product and Manufacturer

Product:	PEXMUX 8 Unit P/N H12592 PEXMUX 16 Unit P/N H12602 PEXMUX Board P/N H12538 CPU-A Board P/N W1522044
Manufacturer:	Bruker BioSpin GmbH Silberstreifen 4 76287 Rheinstetten Germany
Conformity:	EN 61010-1 2010-06 Safety requirements for electrical equipment for measurement, control, and laboratory use - Part 1: General requirements

5.2 General Description

The PEXMUX is available in two configurations:

- PEXMUX 8
- PEXMUX 16

Each PEXMUX unit is built in a 1HE 19" cabinet and includes:

- the PEXMUX board (P/N H12583),
- the CPU-A board (P/N W1522044) and,
- a power supply.

The **PEXMUX 8** Unit (P/N H12592) contains one PEXMUX board and the **PEXMUX 16** Unit (P/N 12602) includes two PEXMUX boards.

The PEXMUX is realized as a splitter to connect one FCTRL to up to 16 SGU's (up to 8 on each PEXMUX board).

Special MRI applications use the PEXMUX unit for parallel excitation sequences.




Figure 5.1: The PEXMUX Unit Wiring

5.3 Features

- The PEXMUX 8 Unit has a 48-Bit LVDS receiver input and 8 x 48-Bit LVDS transmit channels working at 80 MHz.
- The PEXMUX 16 Unit has two 48-Bit LVDS receiver inputs, each with 8 x 48-Bit LVDS transmit channels working at 80 MHz.
- The units have two operation modes:

Parallel Mode: All selected outputs are updated with the input value.

Serial Mode: The selected channels will first be successively preloaded and the updated simultaneously.

- Minimum duration is 25 ns in parallel mode; Maximum duration is 200 ns in sequential mode if 8 channels are in use.
- Minimum resolution is 12.5 ns.
- The time gap between receiving and sending of the first LVDS frame is 225 ns.

5.4 Functional Description

The PEXMUX units receive data from the IPSO FCTRL and distributes the data to the 1–8 LVDS outputs (depending on the configuration). The whole logic is implemented in an FPGA (Altera Cyclone).

The unit is controlled in real time from the pulse program sequence by the register bit (9:0) of the LVDS A-Word and the REG_Val bits (1:0) of the B-Word which are send by the IPSO FCTRL.

Bit	Field	Number	Splitter Information
48	PAR	1	
47	SYNC	1	
46	WID	1	0
45	PLS	3	x
44	PA		x
43	A		x
42	REG 9:0	10	x
41			Update
40...33			8 selection bit
32...17	PHASE 15:00	16	x
16...1	SHAPE 15:0	16	x

Table 5.1 LVDS Word Structure Word A at the LVDS-Interface (Transfer FCTRL to SGU)

Bit	Field	Number	Splitter Information
48	PAR	1	
47	SYNC	1	
46	WID	1	x
45	NCO_SEL 2:0	3	x
44			x
43			x
42	REG_VAL 1:0	2	0
41			0
40	F_VAL 2:0	3	x
39			x
38			x
37	PH_VAL 1:0	2	x
36			x
35	SH_VAL	1	
34...1	F_DATA 33:0	34	x

Table 5.2 LVDS Word Structure Word B at the LVDS-Interface (Transfer FCTRL to SGU)

1. The Register Bit "Reg 9:0" includes the information to preload the selected output register, if Reg_Val 1:0 = 0. If the update bit is set the preloaded register will then be transferred to the LVDS output. The selection is stored until a new mode command with Reg 9:0 ≠ 0 is received. The data transferred via the PEXMUX will have a constant delay of 225 ns independent of the operation mode (parallel or serial).
2. If Reg_Val 1:0 ≠ 0 or Reg_Val 1:0 = 0 and Reg 9:0 = 0 then the received LVDS frame (A-Word, B-Word) is transferred to the preselected LVDS outputs.
3. If the update bit appears before 8 LVDS frames (A-Word, B-Word) are received (serial mode) then the output to the SGU is delayed, so that the time from receive to send is always 225 ns.

The decision table below shows the controlling of the register:

	Reg_Val 1:0	Reg 9	Reg 8	Reg 7:0	Description
		Mode	Update	Select-Bits	
Parallel	≠ 0	x	x	x	In default mode, after power-up, only the register 1 is selected and any incoming data is transferred to Out1. The frame (A/B-Word) is shifted through the register array and transferred to the output after 200 ns. LVDS input (A/B-Word) will be loaded in the preselection register controlled by the command below.
	= 0	x	x	<= 0>	
	= 0	1	x	<selection>	Depending on the <selection> the frame (A/B-Word) is loaded in the appropriate register, parallel shifted through the register array and transferred to the output.
Serial	= 0	0	0	<≠ 0; 1 out n>	The first frame (A/B-Word) is loaded in the first level of the selected register.
	= 0	0	0	<≠ 0; 1 out n>	The second frame is loaded in the first level of the selected register.
	= 0	0	1	<≠ 0; 1 out n>	The third frame is loaded in the first level of the selected register array and then (because the update bit is set) the preloaded words of level 1 are shifted to level 3. The data (all 3 frames) then shifts 5 levels further before sending parallel to the outputs.

Table 5.3 Loading the PEXMUX Register Array

5.5 The PEXMUX Board H12583




Figure 5.2: PEXMUX Board

5.6 JTAG Structure

The board possesses a JTAG interface which can be used for programming and testing.

Connector	St1
JTAG Chains	
Devices	U3: Configuration - EEPROM EPC8
	U6: FPGA EP2C70

Table 5.4 JTAG Structure on the PEXMUX Board

Pin Number	Signal		Pin Number	Signal
1	TRST*		6	JTAG Power +5V
2	Cable Detection		7	TMS
3	TDO		8	GND
4	GND		9	TCK
5	TDI		10	GND




Figure 5.3: JTAG Connector Pin Assignment: Top View Male

Table 5.5 JTAG Connector Pin Assignment

Signal	Description
TRST*	JTAG TAP reset. When asserted low, the TAP controller is asynchronously forced to enter a reset state, which in turn asynchronously initializes other test logic. An unterminated TRST_I produces the same result as if it were driven high. The TAP controller must be reset before the chip can function in normal operating mode.
Cable Detection	Open collector input. Not used in the PEXMUX.
TDO	JTAG serial data out. Signal TDO is the serial output through which test instructions and data from the test logic.
TDI	JTAG serial data in. Signal TDI is the serial input through which JTAG instructions and test data enter the JTAG interface. The new data on TDI is sampled on the rising edge of TCK. An unterminated TDI produces the same result as if TDI were driven high.
TMS	JTAG test mode select. Signal TMS causes state transitions in the test access port (TAP) controller. An un-driven TMS has the same result as if it were driven high.
TCK	JTAG boundary-scan clock. Signal TCK is the clock controlling the JTAG logic.

Table 5.6 IPSO JTAG Signal Description

5.7 Power Requirements

Pin Number	Signal			
1	3.3 V			
2	3.3 V			
3	GND			
4	GND			




Figure 5.4: ST2 Power Connector Pin Assignment: Top View Male

Table 5.7 ST2 Power Connector Pin Assignment

Part Number	Assembly	+3.3 V
PEXMUX 16 Unit H12502	2XPEXMUX Board CPU-A Board	2 Amp

Table 5.8 Part Number for PEXMUX 16 Unit

5.8 Pin Allocation for the LVDS Connectors




Figure 5.5: Pin Assignment for the 48-Bit LVDS Connector at the PCB

5.8.1 CHASSIS

The CHASSIS is a separate plane in the PCB layer stack. This plane is stacked close by the ground plane, giving a very tight capacitive (only capacitive) and low inductance coupling to GND. The chassis plane is screwed together with the external chassis along the front edge near the connectors and the line drivers.

This solution reduces the digital noise at that point and the noise which is picked up by the driver and carried to the outside. In addition this avoids parasitic current through the GND plane which could be caused by potential differences of the remote device.

Function	Wire Type	Transmitter Signal	Receiver Signal	Pin Nr.
Signal: Differential pair of the received serial transmit clock connected to the corresponding inputs of the transmitter.	twisted and shielded	TxCLK_P	RxCLK_P	6
		TxCLK_M	RxCLK_M	18
LVDS Gnd		26		
Shield: Common drain wire of all separate shields, connected to CHASSIS.				
Signal: Differential pair of the received serial data stream connected to the corresponding inputs of the transmitter.	twisted and shielded	TxIN_P0	RxIN_P0	3
		TxIN_M0	RxIN_M0	15
LVDS Gnd		26		
Shield				
Signal	twisted and shielded	TxIN_P1	RxIN_P1	4
		TxIN_M1	RxIN_M1	16
LVDS Gnd		26		
Shield				
Signal	twisted and shielded	TxIN_P2	RxIN_P2	5
		TxIN_M2	RxIN_M2	17
LVDS Gnd		26		
Shield				
Signal	twisted and shielded	TxIN_P3	RxIN_P3	9
		TxIN_M3	RxIN_M3	21
LVDS Gnd		26		
Shield				
Signal	twisted and shielded	TxIN_P4	RxIN_P4	10
		TxIN_M4	RxIN_M4	22
LVDS Gnd		26		
Shield				
Signal	twisted and shielded	TxIN_P5	RxIN_P5	11
		TxIN_M5	RxIN_M5	23
LVDS Gnd		26		
Shield				
Signal	twisted and shielded	TxIN_P6	RxIN_P6	12
		TxIN_M6	RxIN_M6	24
LVDS Gnd		26		
Shield				
Signal	twisted and shielded	TxIN_P7	RxIN_P7	13
		TxIN_M7	RxIN_M7	25
LVDS Gnd		26		
Shield				
USB signal pair, left open.	twisted and shielded	USB+		1
		USB-		14
Shield of the USB signal pair, connected to CHASSIS		USB Gnd		2

Table 5.9 LVDS Connector: Cable and Pin Assignment

Signal: Connected to Bit1 of register „chan-conf“ on F- and G-Controller.	individual	CHANNEL_DETECT0	7
Signal: Connected to Bit0 of register „chan-conf“ on F- and G-Controller.	individual	CHANNEL_DETECT1	20
VCC of USB power, left open.	individual	USB pwr	19
GND of USB power, connected to GND.	individual	USB Gnd	8
Common shield of the entire bundle.	shield	CHASSIS	body

Table 5.9 LVDS Connector: Cable and Pin Assignment



6 Transport, Packaging and Storage

The PEXMUX units contain sensitive and damageable parts and assemblies that must be handled with care.

i Installation and initial commissioning must only be carried out by employees of the manufacturer or persons authorised by the manufacturer.

Nevertheless, it may happen in the course of installation and further use that the owner's operating or maintenance personnel may be entrusted with handling pack units. In this regard, it is vital to observe the instructions listed below.

6.1 Inspection at Delivery

When the PEXMUX is delivered it has to be inspected for completeness and transport damage.

Proceed as follows in the event of externally apparent transport damage:

- Do not accept the delivery, or only accept it subject to reservation.
- Note the extent of the damage on the transport documentation or the shipper's delivery note.
- All the shipping cartons should be stored for further investigation.
- Initiate complaint procedures immediately.

6.2 Packaging

The individual packages are packaged in accordance with anticipated transport conditions. Only environmentally friendly materials have been used in the packaging.

The packaging is intended to protect the individual components from transport damage, corrosion and other damage prior to assembly. Therefore do not destroy the packaging and only remove it shortly before assembly.

6.2.1 Handling Packaging Materials

Dispose of packaging material in accordance with the relevant applicable legal requirements and local regulations.

6.3 Storage

When the PEXMUX is not installed immediately, it has to be stored in the original packing under the following conditions:

- Do not store outdoors.
- Store in dry and dust-free conditions.
- Do not expose to aggressive media.
- Protect against direct sunlight.
- Avoid mechanical shocks.
- Storage temperature: 5°C to 40°C.
- Relative humidity: Maximum of 80% for temperatures up to 31°C, linearly decreasing to 50% at a temperature of 40°C.
- If stored for longer than 3 months, regularly check the general condition of all parts and the packaging. If necessary, top-up or replace preservatives.

6.4 Transport to Manufacturer

When the unit must be returned to the manufacturer for a major repair, use the original packaging for transportation.

Include a good description of the problem.

7 Installation and Initial Commissioning

WARNING



Danger to life from incorrect installation!

Errors during installation may result in potentially hazardous situations and may cause significant damage to property.

- ▶ Installation must only be carried out according to the instructions provided in this manual.

7.1 Safety

Improper initial commissioning

WARNING



Danger of injury from improper initial commissioning!

Improper initial commissioning can result in serious injury and significant damage to property.

- ▶ Before the initial commissioning, ensure that all installation work has been carried out and completed in accordance with the information and instructions in this manual.

Electrical system

WARNING



Danger to life or material damage from electric power!

Contact with live parts may prove fatal. When switched on, electric components are under high voltage.

- ▶ Switch off the power supply before starting work and make sure that it cannot be switched on again.

Securing to prevent restart



WARNING

Danger to life from an un authorised restart!

In the event of an un authorised restart of the power supply during installation, there is a danger of serious injuries or death for persons in the danger zone.

- ▶ Switch off all power supplies before starting work and make sure they cannot be switched on again.

7.2 Installation

All the requirements concerning the environment described in the "[Operating Environment](#)" on page 12 have to be met before the unit can be installed.

To reduce the risk of electric shock and malfunctioning, install these devices in a temperature and humidity controlled indoor area free of conductive contaminants. The power supply cord is intended to serve as the disconnect device. The socket outlet should be near the equipment and easily accessible.

7.2.1 External Connections

The following cables have to be connected before the IPSO can be put into operation:

1. Connect the LAN (Ethernet) RJ45 connector to the Ethernet connector on the HWIDS unit.
2. Connect the LVDS Input (TX-Controller) to an IPSO FCTRL channel. With the PEXMUX 16, connect the second input to a second IPSO FCTRL channel.
3. Connect the LVDS Outputs (up to 8 for the PEXMUX 8, up to 16 for the PEXMUX 16) to the appropriate SGU's.




Figure 7.1: External Connections for the PEXMUX8 Unit

The Ethernet link is used to read the BIS configuration of the unit via the WINBIS and SOAP interface program. The BIS information is stored on the CPU-A board and contains a list of 4 entries:

- PEXMUX Board 1,
- PEXMUX Board 2,
- CPU-A board,
- PEXMUX Unit.

For a detailed description of the operation refer to the PEXMUX_doc_20080813.

To prevent the sending of corrupted data after power up the PEXMUX Unit should be first powered up after the IPSO AQS Unit is ready and connected.

The IPSO Test Program has special commands used to test the correct function of the PEXMUX unit. In this case the LVDS input is connected to the appropriate FCTRL channel and the output is connected to the IPSO RCTRL channel. The PEXMUX outputs are sequentially tested from OUT1 to OUT8.



8 Operating Instructions

The only operations permitted by the user are:

- Connecting the data interface cables (see "Installation" on page 30).
- Switching on the unit using the power switch (see Figure 7.1: for location of the switch).

Once the unit is switched on there are no further operating procedures for the PEXMUX unit. While the unit is running the power switch will be lighted, indicating that the unit is operating.

During operation the user must take great care in following the safety related information provided in this manual.

WARNING



Danger of injury if personnel are insufficiently qualified.

If unqualified personnel perform work on the unit or are in the units danger zone, hazards may arise which can cause serious injury and substantial damage to property:

- ▶ Except for CRU's, all operations inside the PEXMUX must be carried out by a Bruker service engineer or an authorized agent.

8.1 Operator Protection

The electronic circuitry of the PEXMUX operates with low and safe voltages, except for the power supply and its mains connection. Nevertheless, any electrical equipment may become a source of danger under extreme conditions:

WARNING



Risk to life due to extreme weather conditions.

During a lightning storm contact with cables may be life-threatening:

- ▶ Ensure that the electrical power source is properly grounded.
- ▶ Do not loosen, connect or touch any cables during a lightning storm.



WARNING

Material or personnel damage from worn or damaged parts.

Damaged or worn parts may result if one or the following circumstances occurs:

- ▶ The power cord, power plug or power supply are cracked, brittle or damaged.



WARNING

Electrical hazard from electrical shock.

A life threatening shock may result when the housing is open during operation.

- ▶ Disconnect the device from the electrical power supply before opening the device. Use a voltmeter to verify that the device is not under power!
- ▶ Be sure that the power supply cannot be reconnected without notice.
- ▶ The housing must be closed during operation.

8.2 Unit Protection

NOTICE

Material damage hazard from electrostatic discharge (ESD).

Friction between material being conveyed may result in the buildup of electrostatic potential. Contact with the uncovered metal of the printed circuit board (PCB) may result in material damage.

- ▶ Potential equalisation must be ensured before making contact with parts.
- ▶ Where appropriate use ESD flooring and wear ESD shoes.

NOTICE

Material damage hazard from incorrect cable connections.

Incorrect cable connections or disconnecting cables while the unit is under operation may result in material damage.

- ▶ Do not connect a transmitter to the LVDS OUT connectors (OUT1 - OUT8) of the PEXMUX unit!
- ▶ A LVDS cable should never be removed from or connected to an operating controller. Corrupted data may be sampled as valid.

9 Maintenance and Cleaning

9.1 Safety

WARNING



Danger to life or material damage from electric power!

Contact with live parts may prove fatal. When switched on, electric components are under high voltage.

- ▶ Switch off the power supply before starting work and make sure that it cannot be switched on again.

WARNING



Danger to life from an unauthorised restart!

In the event of an unauthorised restart of the power supply during maintenance, there is a danger of serious injuries or death for persons in the danger zone.

- ▶ Switch off all power supplies before starting work and make sure they cannot be switched on again.

WARNING



Danger of injury from improperly executed maintenance work!

Improper maintenance may result in serious injury and significant damage to property.

- ▶ Ensure sufficient assembly space before starting work.
- ▶ Pay attention to orderliness and cleanliness in the assembly location! Loosely stacked or scattered components and tools could cause accidents.
- ▶ If components have been removed, pay attention to correct assembly, refit all fixing elements and comply with bolt tightening torques.

Before the restart, ensure that:

- ▶ All maintenance work has been carried out and completed in accordance with the specifications and instructions in the manual.
- ▶ No persons are in the danger zone.
- ▶ All covers and safety devices are installed and functioning properly.

9.1.1 Environmental Protection

Observe the following environmental protection instructions during maintenance work:

- In respect to all lubrication points supplied manually with lubricant, remove any escaping, used or surplus grease and dispose of it in accordance with applicable local regulations.
- Catch replaced oils in suitable containers and dispose of in accordance with applicable local regulations.

9.2 Maintenance

The PEXMUX unit can be expected to have a long and trouble-free life with a minimum of preventive maintenance. Environmental issues are essential in determining the reliability. The temperature and humidity have to be within specifications. All servicing must be performed by qualified service personnel.

Before maintenance, repair or shipment, the unit must be completely switched off and unplugged or disconnected and dismantled from its rack.

9.3 Cleaning

Cleaning the surface of the enclosure and/or front panel can be carried out by the customer, if the following instructions are adhered to.

1. Switch off the equipment and unplug the power cable and all data cables.
2. Clean up the surface with a dry or damp cloth.
3. Let the enclosure completely dry before installing
4. Connect all cables and power up.

9.4 Service Requests

In case of questions or problems, please contact your nearest Bruker office or representative. See ["Contact" on page 39](#) for details.

10 Dismantling and Disposal

After the lifespan of the instrument, Bruker takes responsibility for disassembly and disposal in accordance with the European directive 2012/19/EC WEEE. Bruker BioSpin GmbH offers to take back the components free of charge after usage at the customer site upon request by the customer. If the customer wants to arrange disposal on their own, then this has also to be stated when the product is ordered.



11 Contact

Manufacturer:

Bruker BioSpin NMR

Silberstreifen

D-76287 Rheinstetten

Germany

Phone: +49 721-5171-02

<http://www.bruker.com>

WEEE DE43181702

NMR Hotlines

Contact our NMR service centers.

Bruker BioSpin NMR provide dedicated hotlines and service centers, so that our specialists can respond as quickly as possible to all your service requests, applications questions, software or technical needs.

Please select the NMR service center or hotline you wish to contact from our list available at:

<http://www.bruker.com/service/information-communication/helpdesk.html>



Appendix A

A.1 Warning Signs

Danger

Key Word and Symbol	6
---------------------------	---

Warning

Danger of injury from improper initial commissioning!	29
Danger of injury from improperly executed maintenance work!	35
Danger of injury if personnel are insufficiently qualified.	33
Danger to life from an unauthorised restart!	35
Danger to life from incorrect installation!	29
Danger to life or material damage from electric power!	29, 35
Electrical hazard from electrical shock.	34
Key Word and Symbol	6
Material or personnel damage from worn or damaged parts.	10, 34
Risk to life due to extreme weather conditions.	33

Caution

Key Word and Symbol	6
---------------------------	---

Note

Installation and initial commissioning must only be carried out by... ..	27
Only trained service personnel are allowed to mount... ..	10

Notice

Key Word and Symbol	6
Material damage hazard from electrostatic discharge (ESD).	34
Material damage hazard from incorrect cable connections.	34

Safety Instructions

Key Word and Symbol	6
---------------------------	---



A.2 Figures

Figure 5.1:	The PEXMUX Unit Wiring	16
Figure 5.2:	PEXMUX Board	20
Figure 5.3:	JTAG Connector Pin Assignment: Top View Male	21
Figure 5.4:	ST2 Power Connector Pin Assignment: Top View Male.....	22
Figure 5.5:	Pin Assignment for the 48-Bit LVDS Connector at the PCB	23
Figure 7.1:	External Connections for the PEXMUX8 Unit	31



A.3 Tables

Table 5.1	LVDS Word Structure Word A at the LVDS-Interface (Transfer FCTRL to SGU) 17	
Table 5.2	LVDS Word Structure Word B at the LVDS-Interface (Transfer FCTRL to SGU) 18	
Table 5.3	Loading the PEXMUX Register Array	19
Table 5.4	JTAG Structure on the PEXMUX Board	20
Table 5.5	JTAG Connector Pin Assignment.....	21
Table 5.6	IPSO JTAG Signal Description.....	21
Table 5.7	ST2 Power Connector Pin Assignment.....	22
Table 5.8	Part Number for PEXMUX 16 Unit	22
Table 5.9	LVDS Connector: Cable and Pin Assignment.....	24



Revision History List

Index	Date	Alteration Type



Bruker Corporation

info@bruker.com
www.bruker.com

Order No: Z31929